


## Advanced Wound Care



### Antibacterial **MEDIHONEY**<sup>®</sup> *Apinate*<sup>™</sup> Dressing

MediHoney<sup>®</sup> Antibacterial Honey Apinate<sup>™</sup> dressing combines Medical Grade Manuka Honey with a calcium alginate.

**95%**  
ANTIBACTERIAL  
MANUKA HONEY\*

#### Features and Benefits

- Made with Medical Grade Manuka Honey
- The calcium alginate provides additional absorbency capability
- Delivers broad-spectrum antimicrobial activity<sup>1</sup>
- MediHoney aids and supports autolytic debridement and a moist wound healing environment in acute and chronic wounds and burns<sup>2-5</sup>
- MediHoney has a low pH of 3.5-4.5
- Lowering pH has been associated with wound healing benefits
- Helps to reduce malodour<sup>6</sup>
- Has been safely used on patients of all ages
- Can be easily cut, shaped or moulded to fit any wound type including cavities when rope version is used
- Dressing change frequency will depend on the condition of the patient as well as the level of wound exudate

# MEDIHONEY® Apinate™

## Antibacterial Dressing

### MediHoney® Apinate™ Antibacterial Dressing

	Reference	Description	Units/Box	PIP Code	NHS Code
Rope	793	1.9 cm x 30 cm	5/box	366-5429	EJE077
	794	5 cm x 5 cm	10/box	366-5411	EJE076
Dressing	795	10 cm x 10 cm	5/box	347-0929	EJE000



Dressing




Rope



Leg Ulcer



Pressure Ulcer



Infected Wound

#### Indications

- Leg ulcers (venous, arterial and mixed aetiology ulcers)
- Diabetic foot ulcers
- Pressure ulcers
- Infected and malodorous wounds
- Donor sites and recipient graft sites
- Burns
- Surgical wounds

#### Contraindications

- Third-degree burns
- Patients with a known sensitivity to alginates or honey
- To control heavy bleeding

**References:** 1. Maddocks S, Lopez M, Rowlands R, Cooper R. Manuka honey inhibits the development of *Streptococcus pyogenes* biofilms and causes reduced expression of two fibronectin-binding proteins. *Microbiology*. 2012;781-790. 2. Regulski, M. A novel wound care dressing for chronic leg ulcerations. *Podiatry Management*. 2008; November: 235-246. p. 235-246. 3. Robson, V., Dodd, S and Thomas, S. Standardized antibacterial honey (MediHoney®) with standard therapy in wound care: randomized clinical trial. *Journal of Advanced Nursing*, 2009;565-575. 4. Bateman S, Graham T. The Use of MediHoney® Wound Gel on surgical wounds post-CABG. *WOUNDS UK*; 2007; (3);76-83. 5. Cadogan, J. The use of honey to treat an ulcer on the heel of a person with diabetes. *The Diabetic Foot Journal*;2008; (1);43-45. 6. Robson V, Yorke J, Sen R, Lowe D, Rogers S. Randomised controlled feasibility trial on the use of medical grade honey following microvascular free tissue transfer to reduce the incidence of wound infection. *The British J Oral & Maxillofacial Surgery*. 2012;50(4):321-327.

**Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.**

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

#### Additional information for EMEA Customers only:

Products mentioned in this document are CE class I, IIa, IIb or III devices. Contact Integra should you need any additional information on devices classification. All the medical devices mentioned on this document are CE marked according to European council directive 93/42/EEC on medical devices and its relatives, unless specifically identified as "NOT CE MARKED".

### For more information or to place an order, please contact:

#### Customer Service

United Kingdom +44 (0) 1264 312 725 ■ +44 (0) 1264 312 821 fax  
[integrallife.eu](http://integrallife.eu)

#### Manufacturer:



Derma Sciences, Inc.  
 104 Shorting Road  
 Toronto, Ontario M1S 3S4 ■ Canada




Integra LifeSciences Services (France)  
 Immeuble Séquoia 2  
 97 Allée Alexandre Borodine  
 Parc Technologique de la Porte des Alpes  
 69800 Saint Priest – France



**MEDIHONEY®**  
Antibacterial Dressing

## Advanced Wound Care



### Antibacterial **MEDIHONEY®** *Apinate™ Dressing*

MediHoney® Antibacterial Honey Apinate™ dressing combines Medical Grade Manuka Honey with a calcium alginate.

**95%**  
ANTIBACTERIAL  
MANUKA HONEY\*

#### Features and Benefits

- Made with Medical Grade Manuka Honey
- The calcium alginate provides additional absorbency capability
- Delivers broad-spectrum antimicrobial activity<sup>1</sup>
- MediHoney aids and supports autolytic debridement and a moist wound healing environment in acute and chronic wounds and burns<sup>2-5</sup>
- MediHoney has a low pH of 3.5-4.5
- Lowering pH has been associated with wound healing benefits
- Helps to reduce malodour<sup>6</sup>
- Has been safely used on patients of all ages
- Can be easily cut, shaped or moulded to fit any wound type including cavities when rope version is used
- Dressing change frequency will depend on the condition of the patient as well as the level of wound exudate

**INTEGRA®**  
LIMIT UNCERTAINTY

\* by weight

## MediHoney® Apinate™ Antibacterial Dressing

	Reference	Description	Units/Box	PIP Code	NHS Code
<b>Rope</b>	793	1.9 cm x 30 cm	5/box	366-5429	EJEO77
	794	5 cm x 5 cm	10/box	366-5411	EJEO76
<b>Dressing</b>	795	10 cm x 10 cm	5/box	347-0929	EJEO00



Dressing




Rope



Leg Ulcer



Pressure Ulcer



Infected Wound

### Indications

- Leg ulcers (venous, arterial and mixed aetiology ulcers)
- Diabetic foot ulcers
- Pressure ulcers
- Infected and malodorous wounds
- Donor sites and recipient graft sites
- Burns
- Surgical wounds

### Contraindications

- Third-degree burns
- Patients with a known sensitivity to alginates or honey
- To control heavy bleeding

**References:** 1. Maddocks S, Lopez M, Rowlands R, Cooper R. Manuka honey inhibits the development of *Streptococcus pyogenes* biofilms and causes reduced expression of two fibronectin-binding proteins. *Microbiology*. 2012;781-790. 2. Regulski, M. A novel wound care dressing for chronic leg ulcerations. *Podiatry Management*. 2008; November: 235-246. p. 235-246. 3. Robson, V., Dodd, S and Thomas, S. Standardized antibacterial honey (MediHoney®) with standard therapy in wound care: randomized clinical trial. *Journal of Advanced Nursing*, 2009;565-575. 4. Bateman S, Graham T. The Use of MediHoney® Wound Gel on surgical wounds post-CABG. *WOUNDS UK*; 2007; (3);76-83. 5. Cadogan, J. The use of honey to treat an ulcer on the heel of a person with diabetes. *The Diabetic Foot Journal*;2008; (1);43-45. 6. Robson V, Yorke J, Sen R, Lowe D, Rogers S. Randomised controlled feasibility trial on the use of medical grade honey following microvascular free tissue transfer to reduce the incidence of wound infection. *The British J Oral & Maxillofacial Surgery*. 2012;50(4):321-327.

**Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.**

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

### Additional information for EMEA Customers only:

Products mentioned in this document are CE class I, IIa, IIb or III devices. Contact Integra should you need any additional information on devices classification. All the medical devices mentioned on this document are CE marked according to European council directive 93/42/EEC on medical devices and its relatives, unless specifically identified as "NOT CE MARKED".

For more information or to place an order, please contact:

### Customer Service

United Kingdom +44 (0) 1264 312 725 ▪ +44 (0) 1264 312 821 fax  
[integralife.eu](http://integralife.eu)

### Manufacturer:



Derma Sciences, Inc.  
 104 Shorting Road  
 Toronto, Ontario M1S 3S4 ▪ Canada



Integra LifeSciences Services (France)  
 Immeuble Séquoia 2  
 97 Allée Alexandre Borodine  
 Parc Technologique de la Porte des Alpes  
 69800 Saint Priest – France